



HEALTH CARE PRODUCTS



In a Heartbeat

Good Health, Good Choices, Good Living



TABLE OF CONTENTS

• Pressure machines	4-7
Digital blood pressure monitor (arm type)	4-5
Aneroid sphygmomanometers	6-7
Child aneroid sphygmomanometer	7
• Spare parts accessories	8-9
• Blood glucose meter	10
• Nebulizers	11-14
Compressor nebulizer	11-14
Dolphin compressor nebulizer for kids	14
• Spacers for aerosol	14
• Oxygen concentrator (5 Liters + 10 Liters)	15
• Oximeters	16
• Electronic scales	17-19
Electronic personal scales	17
Body composition scales	18
Height & weight baby scale	19
• Thermometers	20-22
Digital thermometers	20
Ear & Infrared thermometers	21-23
• Alcohol, Hydrogene peroxyde solution	24-25
• Hygienic masks and hand sanitizers	26-33
• Rubber hot water bag	34
• Gauze, syringes, bandages and First aid kit	35-57
• Cotton	58-60
• Ovulation and pregnancy tests	61
• Adults pads	62
• Wipes	63-64
• Skin care products	65
• Sunscreen products	66

For A Healthy Heart Rate



DIGITAL BLOOD PRESSURE MONITORS

Arm Type Automatic



Item #	Description	Qty/Pack	Barcode
OPM-1317	Digital Blood Pressure Monitor	1Pc	 7 261894 213171



Item #	Description	Qty/Pack	Barcode
OPM-1339	Digital Blood Pressure Monitor	1Pc	 7 261894 213393

Digital BPM	OPM - 1317	OPM - 1339	OPM - 1356
LCD display size	Large	Medium	Medium
Date/time display	✓	✓	✓
Low battery indication	✓	✓	✓
Cuff check indicator	X	X	✓
Movement detection	X	X	✓
Number of users	2	2	2
Memory set/user	90 sets	90 sets	90 sets
SBP, DBP & pulse rate display	✓	✓	✓
Alarm clock reminding function	X	X	X
USB cable	✓	✓	✓
Stand	X	X	X
Case for cuff	X	✓	✓
IHB detector	✓	✓	✓
Inbuilt battery	X	X	X

- Maximum accuracy
- Irregular heart beat detection
- 85-90% less power consumption



Item #	Description	Qty/Pack	Barcode
OPM-1356	Digital Blood Pressure Monitor	1Pc	 7 261894 213560

DIGITAL BLOOD PRESSURE MONITOR

Arm Type Automatic

Digital BPM	OPM - 1363	OPM - 1310
LCD display size	Medium	Large
Date/time display	✓	✓
Low battery indication	✓	✓
Cuff check indicator	✓	X
Movement detection	✓	X
Number of users	2	2
Memory set/user	90 sets	60 sets
SBP, DBP & pulse rate display	✓	✓
Alarm clock reminding function	X	✓
USB cable	✓	X
Stand	X	✓
Case for cuff	✓	X
IHB detector	✓	✓
Inbuilt battery	X	✓

Item #	Description	Qty/Pack	Barcode
OPM-1363	Digital Blood Pressure Monitor	1Pc	 7 261894 213638
OPM-1310	Digital Blood Pressure Monitor	1Pc	 7 261894 213102



OPM - 1363



OPM - 1310

ANEROID SPHYGMOMANOMETER

With Black Nylon Cuff

- Measuring range: 0-300mmHg
- Square shaped dial 14mm x 145mm gradient adjustable for easy reading
- Standard cotton cuff
- 2-tube adult latex blade
- Standard latex bulb
- Air-release valve with spring
- Standard end valve
- Spiral PVC tube (1.8m)

Measuring Method	Oscillometry
Measuring Range	0-300mmHg (blood pressure)
Precision	±3mmHg
Scale Graduation	2mmHg



- Manually-operated
- Measurements are taken from the upper arm

Item #	Description	Qty/Pack	Barcode
OPM-1324	Aneroid Sphygmomanometer	1Pc	 7 261894 213249

ANEROID SPHYGMOMANOMETER

With Black Nylon Cuff



CHILD ANEROID SPHYGMOMANOMETER

Latex Free



Item #	Description	Qty/Pack	Barcode
OPM-1311	Aneroide Sphygmomanometer	1Pc	 726189421311

Item #	Description	Qty/Pack	Barcode
OPM-1314	Child Aneroide Sphygmomanometer	1Pc	 726189421314

SPARE PARTS ACCESSORIES

NEBULIZING KIT

Disposable nebulizing kit consists of:

Silicone mask | Mouthpiece | Nebulizer chamber | Tubing | Elastic band



Small size



Medium size



Large size

Directions for use:

1. Open the package within the validity, take out the disposable nebulizing kit.
2. Unscrew the lid of the nebulizer chamber and place the prescribed medication.
3. Reattach the lid and close the nebulizer chamber.
4. Put the mask on the patient's face, centering both nose and mouth.
5. Gently adjust the ends of the elastic strap to secure the mask.
6. Connect nebulizing kit to the compressor, start the inhalation.
7. Adjust the mist volume by turning the knob right or left over the nebulizer chamber cap.

Item #	Description	Qty/Pack	Barcode
OPM-1346	Nebulizing kit small size	1 Pc	 7 261894 213461
OPM-1347	Nebulizing kit medium size	1 Pc	 7 261894 213478
OPM-1348	Nebulizing kit large size	1 Pc	 7 261894 213485

NEBULIZER MASK SET FOR:

- Adults
- kids



DISPOSABLE OXYGEN MASK



Item #	Description	Qty/Pack	Barcode
OPM-1328	Nebulizer mask set for adult	1Pc	7 261894 213287
OPM-1326	Nebulizer mask set for Kids	1Pc	7 261894 213263
OPM-1345	Disposable Oxygen Mask	1Pc	7 261894 213454

REPLACEMENT BLOOD PRESSURE CUFF



- Complete home blood pressure kit allows for easy monitoring by an individual or with assistance from another person.
- Small adult (22-32cm).

Item #	Description	Qty/Pack	Barcode
OPM-1308	Replacement Blood Pressure Cuff	1Pc	7 261894 213089

BLOOD PRESSURE CUFF AND TWO-TUBE BLADDER

Item #	Description	Qty/Pack	Barcode
OPM-1309	Blood Pressure Cuff and Two-Tube Bladder	1Pc	7 261894 213096

- Double tube bladders.
- Blue nylon cuff.
- Adult size.



BLOOD GLUCOSE METER

- Gold AQ Plus strips stable in test
- Fast and accurate
- Testing result in 5 seconds
- Virtually pain-free
- For capillary blood and venous whole blood

Attention:

- Only use with Optimal® Gold AQ Plus Blood glucose test strip
- For in vitro diagnostic use only

0.8 Micro Liters



5 Seconds

NO CODING

BLOOD GLUCOSE TEST STRIPS

99.9% GOLD ELECTRODE



GOLD+

Wide HCT Range



BLOOD LANCETS



- Sterile
- Ultra fine, tri-bevel tip for virtually painless sampling
- Consistent depth penetration
- Universal design fits almost all lancing devices
- 28 Gauge

Warning:

- Do not use the lancet if the protective cap has been previously removed.
- Do not reuse lancets.
- Do not store lancets in lancing devices.

Item #	Description	Qty/Pack	Barcode
OPM-1352	Blood Glucose Meter	1 pc	7 ¹ 261894 ¹ 213522
OPM-1353	Blood Glucose Test Strips	50 pcs	7 ¹ 261894 ¹ 213539
OPM-1355	Blood Lancets	50 pcs	7 ¹ 261894 ¹ 213553

COMPRESSOR NEBULIZERS

Low Noise, Light Weight

COMPRESSOR NEBULIZERS	OPM - 1349	OPM - 1323
Size	162mm x 132mm x 99mm	150mm x 165mm x 110mm
Rated Voltage	AC 230V, 50HZ	AC 230V, 50HZ
Nebulizing Pressure	60 Kpa ~ 180 Kpa	100 Kpa ~ 190 Kpa
Nebulizing Flow	4 L/min	4 L/min ~ 7L/min
Extreme Pressure	210 Kpa ~ 400 Kpa	>205 Kpa
Rated Power	190 W	60 W
Capacity Of Medication Tank	10 ml	10 ml



OPM - 1349



OPM - 1323



Coughing,
especially at night



Wheezing



Shortness
of breath



Chest tightness,
pain, or pressure

Item #	Description	Qty/Pack	Barcode
OPM-1349	Compressor Nebulizer	1Pc	 7 261894 213492
OPM-1323	Compressor Nebulizer	1Pc	 7 261894 213232

COMPRESSOR NEBULIZERS

For Your Family Comfort

COMPRESSOR NEBULIZERS	OPM - 1321	OPM - 1318
Size	280mm x 180mm x 88mm	330mm x 200mm x 138mm
Rated Voltage	AC 230V, 50HZ	AC 230V, 50HZ
Nebulizing Pressure	60 Kpa ~ 180 Kpa	100Kpa ~ 190Kpa
Nebulizing Flow	≥7L/min	4 L/min ~ 7L/min
Extreme Pressure	210 Kpa ~ 400 Kpa	>205 Kpa
Rated Power	132 W	60 W
Capacity Of Medication Tank	10 ml	10 ml



OPM - 1321



Coughing,
especially at night



Wheezing



Shortness
of breath



Chest tightness,
pain, or pressure

Item #	Description	Qty/Pack	Barcode
OPM-1321	Compressor Nebulizer	1Pc	 7 ² 261894 ² 213218
OPM-1318	Compressor Nebulizer	1Pc	 726189421318



OPM - 1318

WEARABLE MESH NEBULIZER

Small and light | hands free

COMPRESSOR NEBULIZERS	OPM - 1357
Size	54mm x 36mm x 44mm
Power Supply	DC 5 V, 0.5 A
Noise Level	<50 dB
Cleaning Mode	✓
Fully Charged Battery Life	>1 hour
Rated Power	≤3W
Capacity Of Medication Tank	6 ml

- Cleaning Mode
- Ultrasonic Frequency: Approx. 110 khz
- Power Supply: Built-in Lithium Battery (3.7V)
- Rated Power: ≤3W
- Mmad: <5Mm - Fpf: ≥60%
- Nebulizing Rate: ≥0.2 ml/min
- Nebulizer Head Life: 180 Hours
- Fully Charged Battery Life: >1 Hour
- Max./ Min. Liquid Volume: 6 ml/2 ml
- Automatic Shutdown Without Medicine
- Main Unit Size: 54x36x44 Mm | Weight: Approx. 45g



Accessories are available



Item #	Description	Qty/Pack	Barcode
OPM-1357	Wearable Mesh Nebulizer	1Pc	 7 261894 213577

DOLPHIN COMPRESSOR NEBULIZER FOR KIDS

For kids, Low noise, Light weight



- Concise and fluent design
 - Nebulizing flow 4~8 L/min
 - One-button operation for easy use
 - Highly effective delivery of medicine
 - Convenient to clean the dust filter
 - Extreme pressure more than 196Kpa
 - Small enough for home and clinical use
 - Clever handle design to be compatible with medical cup
-
- Deep penetration of nasal and chest congestion
 - Simple 5 - 15 minutes treatment
 - Doctor recommended

Item #	Description	Qty/Pack	Barcode
OPM-1322	Dolphin Compressor Nebulizer	1Pc	

SPACER FOR AEROSOL - 175 ml

For Single Patient Use Only

- Anti-static spacer
- 100% Latex free
- With a facemask for babies
- Creates a secure seal that helps parents administer aerosol medications to infants
- Helps deliver medicine straight to the lungs
- Easily disassembled and reassembled for cleaning



Item #	Description	Qty/ml	Barcode
OPM-1340	0-18 months	175 ml	
OPM-1341	1-5 years	175 ml	
OPM-1342	5 years +	175 ml	

OXYGEN CONCENTRATOR

Oxygen Concentrator	OPM - 1350	OPM - 1351
Mode of use	Single	Dual
Capacity	5 Liters	10 Liters
High oxygen concentration (93%)	✓	✓
Oxygen Purity: $\geq 93\% \pm 3(V/V)$	✓	✓
Stable rotary flowmeter	✓	✓
Home & hospital use	✓	✓
CE/ISO approved	✓	✓
Input power	300W	550W
Voltage	220V/50Hz	220V/50Hz
Flow rate	0-5L/min	0-10L/min
Output pressure	20Kpa-60Kpa	20Kpa-60Kpa
Safety class	Class II, Type A	Class II, Type A



Oxygen purity:
93%±3%



Suitable for
all ages



Fine atomized
particles



OPM - 1350



OPM - 1351

WARNING:

- Do not use near heat or open flames.
- Electric shock hazard, do not remove cabinet.
- Refer all services to qualified service personnel.
- Do not smoke near patient or device.

Item #	Description	Qty/Pack	Barcode
OPM-1350	Oxygen Concentrator (5Liters)	1Pc	 7 ²⁶¹⁸⁹⁴ 213508

Item #	Description	Qty/Pack	Barcode
OPM-1351	Dual Oxygen Concentrator (10Liters)	1Pc	 7 ²⁶¹⁸⁹⁴ 213515

FINGERTIP PULSE OXIMETER

Quick and accurate

Optimal Fingertip Pulse Oximeter is a compact, simple, reliable, and durable monitoring device suitable in a broad range of clinical settings including general practice & home use for all family members.

- Oxygen level
- Pulse rate
- High accuracy
- Lightweight and easy to carry
- Automatic switch on/off
- Intelligent power saving design
- Anti-motion technology
- White LED Light



High performance:

Optical detection technology and anti-slip silicone mat finishing inside enables you more accurate measurement.

Item #	Description	Qty/Pack	Barcode
OPM-1338	Fingertip Pulse Oximeter	1Pc	 7 261894 213386

CHILD FINGERTIP PULSE OXIMETER

Quick and accurate

This small, non-invasive device is designed to measure the amount of oxygen in your child's blood and monitor their heart rate. Using a pulse oximeter is easy and painless—it simply clips onto a fingertip or toe and gives you quick readings.

- One button operation
- Lightweight and portable
- Automatic shut down
- Cute appealing design
- Alarm function
- 20 sets of memory
- Alarm function
- Anti-slip soft silicone
- 5 levels of brightness adjustment
- Multi-direction display
- Test DEMO



Get accurate results to help keep your little one healthy and safe.

Item #	Description	Qty/Pack	Barcode
OPM-1365	Child Fingertip Pulse Oximeter	1Pc	 7 261894 213652

ELECTRONIC PERSONAL SCALES

ULTRA SLIM - 180 Kg capacity

ELECTRONIC SCALES	OPM - 1301 OPM - 1325	OPM - 1334 OPM - 1335
Type	Stainless Steel	Glass
Capacity	180 kg	180 kg
LCD display	✓	✓
Division	100g / 0.2lb	100g / 0.2lb
Low battery, overload indication	✓	✓
Power	2 pcs 1.5V AAA (included)	2 pcs 1.5V AAA (included)
Auto on	✓	✓



OPM - 1301

OPM - 1325



OPM - 1334

OPM - 1335

Item #	Description	Qty/Pack	Barcode
OPM-1301	Bronze Stainless steel Electronic scale	1Pc	 7 261894 213010
OPM-1325	Silver Stainless steel Electronic scale	1Pc	 7 261894 213256

Item #	Description	Qty/Pack	Barcode
OPM-1334	Black Electronic Personal scale	1Pc	 7 261894 213348
OPM-1335	Magenta Electronic Personal scale	1Pc	 7 261894 213355

BODY COMPOSITION SCALES

ULTRA SLIM - 150 Kg capacity

OPM - 1360	OPM - 1361	OPM - 1362
Body composition scale	Body fat and hydration scale	Body fat and hydration scale using BIA principle
Fat percentage	Fat percentage	Fat percentage
Water percentage	Water percentage	Water percentage
Muscles, bones percentage	X	X
BMI feature	BMI feature	BMI feature
Touch buttons	Non-tactile buttons	Non-tactile buttons
Auto user Recognition	Auto user Recognition	Auto user Recognition
12 memories	12 memories	12 memories
White backlight LCD	X	X



OPM - 1360



OPM - 1362



OPM - 1361



OPM - 1361

BMI Body mass index

Auto-Recognition:
Recognize user code according to his weight

Item #	Description	Qty/Pack	Barcode
OPM-1360	Body composition scale	1Pc	 7 261894 213607
OPM-1362	Body Fat & Hydration Scale using BIA	1Pc	 7 261894 213621

Item #	Description	Qty/Pack	Barcode
OPM-1361	Body Fat And Hydration Scale	1Pc	 7 261894 213614

HEIGHT AND WEIGHT BABY SCALE

25 Kg capacity

ALL-IN-ONE
DIGITAL
BABY SCALE

- Ergonomic cradle designed platform
- Infant height measurement with extensible/adjustable platform length
- Maximum 8 users supported, history data review
- Measurements comparison versus previous results
- Easy to read extra large display for height & weight
- Zero button subtracts the weight of diapers or clothing
- Equipped with auto-off technology to conserve battery life



<p>Max. 25kg 55lb</p>	<p>Weight Division 5g/0.01lb</p>	<p>LCD White backlight</p>	<p>LCD SIZE 75x30mm</p>
<p>4PCS (excluded)</p>	<p>8 Memory</p>	<p>Height 46-80cm/ 18.1-31.5in</p>	<p>Height Division 1mm/0.1in</p>

Item #	Description	Qty/Pack	Barcode
OPM-1329	Height & Weight Digital Baby Scale	1Pc	<p>7 261894 213294</p>



DIGITAL THERMOMETERS

OPM - 1332	OPM - 1331
Fixed Tip	Flexible Tip
Result in 30 sec.	Result in 30 sec.
With memory function	With memory function
With beeper	With beeper
Auto power off	Auto power off

OPM - 1333	OPM - 1330
Flexible Tip	Flexible Tip
Result in 30 sec.	Result in 60 sec.
With memory function	With memory function
With beeper	With beeper
Auto power off	Auto power off



Item #	Description	Barcode
OPM-1331	Digital flexible Thermometer	7 261894 213317
OPM-1332	Digital Thermometer	7 261894 213324
OPM-1333	Digital flexible Thermometer	7 261894 213331
OPM-1330	Digital flexible Thermometer	7 261894 213300



OPM - 1332

Fixed tip



OPM - 1331

Flexible tip



OPM - 1333

Flexible tip



OPM - 1330

Flexible tip


PREHEAT Ear Thermometer

- 360° flexible probe
- Preheating sensor technology
- Night mode: Silent temperature measurement with night light
- Age mode: Can be used by babies, children and adults
- Stores 32 sets of memories
- 20 individually packed probe covers
- Probe cover automatically ejected

SPECIFICATIONS	OPB - 1366
Measurement	Body temperature
Mode	Contact type
Alert	Fever alert colors
Memory	32 measurements
LCD	3 colors
Unit	Celsius & Fahrenheit
Automatic power off	60 seconds
Result	in 1 second



OPM - 1366



THE PREHEAT SENSOR TECHNOLOGY
LOWER THE COOLING EFFECT AND
PROMOTES ACCURACY



OPM - 1367



Flexible probe



Night mode

CONTAINS

- Ear Thermometer 1pc
- Thermometer Case 1pc
- Instruction Manual 1pc
- Quick Reference Guide 1pc
- AA alkaline Battery 2pcs
- Probe Cover 20pcs (1 box)



Prevents cross infection

x 20
Probe Covers

Super thin for accurate readings

EAR, FOREHEAD Infrared Thermometer

- Accurate results
- 1 second testing
- Easy to use with 1 single button
- 2 operating modes: Ear and forehead
- Stores the last recorded result
- Three-color backlight to indicate the temperature range
- Ability to switch between Celsius and Fahrenheit display
- Automatic shut-off after 35 seconds to conserve battery life
- Practical, easy to use

SPECIFICATIONS	OPB - 1364
Measurement	Body temperature
Mode	Infrared detection system
Alert	Fever alert colors
Memory	1 last measurement
LCD	3 colors
Unit	Celsius & Fahrenheit
Automatic power off	35 seconds
Result	in 1 second



Ear mode



Forehead mode



Green

Normal Temperature



Yellow

Mild Fever



Red

High Fever



OPM - 1364

DESIGNED FOR THE WHOLE FAMILY

NON-CONTACT INFRARED THERMOMETERS 6 IN 1

Body, Surface, Ambient Temperature, Alarm, Memories 32 Test Data, Three Colors Backlit



**ONE SECOND
MEASURING**

Non-contact Infrared Thermometer	OPM - 1336	OPM - 1327
6 in 1	✓	✓
3~5cm (1.2~2 In) distance from forehead	✓	✓
Audible alarm if temperature is more than 38°C (100.4°F)	✓	✓
Memorize the last 32 temperature measurements	✓	✓
Three colors backlit LCD	✓	✓
Temperature unit can be displayed in either celsius or fahrenheit	✓	✓
Automatic power off	≤30 seconds	≤30 seconds
Longevity use	100.000 Readings	40.000 Readings



OPM - 1336



OPM - 1327

Aim the thermometer probe at the center of the forehead, keep the distance between 3cm and 5cm, and do not touch the forehead directly.

Item #	Description	Qty/Pack	Barcode
OPM-1336	Non-contact infrared thermometer	1Pc	 7 261894 213362
OPM-1327	Non-contact infrared thermometer	1Pc	 7 261894 213270



Quiet mode for sleeping babies



Non-contact measuring



10°C ~ 40°C
(50°F ~ 104°F)



Wipe any sweat



No gear

ALCOHOL: CLEANS AND PROTECTS

70% u.s.p. by volume



OPH - 026

OPH - 051

OPH - 021

OPH - 019

OPH - 017

OPH - 024

OPH - 018

Item #	Description	Qty/ml	Barcode
OPH-026	Alcohol 70%	125ml	
OPH-051	Alcohol 70%	150ml	
OPH-021	Alcohol 70%	200ml	
OPH-019	Alcohol 70%	250ml	

Item #	Description	Qty/ml	Barcode
OPH-017	Alcohol 70%	250ml	
OPH-024	Alcohol 70%	500ml	
OPH-018	Alcohol 70%	1000ml	

75% u.s.p. by volume



OPH - 050

OPH - 052

95% u.s.p. by volume



OPH - 031

OPH - 048

OPH - 020

HYDROGENE PEROXYDE SOLUTION

10 VOL.



OPH - 049

Item #	Description	Qty/ml	Barcode
OPH-050	Alcohol 75%	100ml	7 261894 300505
OPH-052	Alcohol 75%	500ml	7 261894 300529

OPH-031	Alcohol 95%	250ml	7 261894 300031
OPH-048	Alcohol 95%	500ml	7 261894 300048
OPH-020	Alcohol 95%	1000ml	7 261894 300208
OPH-049	Hydrogene Peroxyde Solution	150ml	7 261894 300499

MEDICAL DISPOSABLE FACE MASK

3 Layers 25+25+25GSM

- Three layers design
- Anti dust & fluid
- Comfortable & breathable
- Adjustable nose bridge
- Without glass fibers
- Hypoallergenic
- Very low resistance to breathing
- Nose bar adaptable
- High filtration capacity
- Perfect fitting



Item #	Description	Qty/Pack	Barcode
OPH-371	Medical Disposable Face Mask	50Pcs	 7 261894 303711

MEDICAL PROTECTIVE FACE MASK | KN95

5 layers filtering method for maximum protection

- Disposable
- Nose tip shaping
- Fit a variety of face shape
- 3D stereo curved design
- Smooth breathing
- Skin-friendly fabric
- 99% Soft and comfortable



Item #	Description	Qty/Pack	Barcode
OPH-413	Medical Protective Face Mask	10Pcs	 7 261894 304138

Item #	Description	Qty/Pack	Barcode
OPH-419	Medical Protective Face Mask	50Pcs	 7 261894 304190

KIDS DISPOSABLE FACE MASK

3 protective Layers

- Layer 1: non woven | layer 2: meltblown | layer 3: non woven
- Sized for kids, making it comfy and fun!
- Light breathable material
- Comfortable light earloop

Item #	Description	Qty/Pack	Barcode
OPH-416	Kids Disposable Face Mask	10Pcs	 7 261894 304169



KIDS 3D DISPOSABLE FACE MASK

- Small size: 1 - 4 years | 4 protective layers

- Layer 1 & 2: non woven | layer 3: meltblown | layer 4: non woven
- Sized for kids, making it comfy and fun!
- Light breathable material
- Comfortable light earloop



Item #	Description	Qty/Pack	Barcode
OPH-414	Kids 3D Disposable Face Mask (small)	20Pcs	 7 261894 304145

KIDS 3D DISPOSABLE FACE MASK

- Small size: 4 years + | 4 protective layers

- Layer 1 & 2: non woven | layer 3: meltblown | layer 4: non woven
- Sized for kids, making it comfy and fun!
- Light breathable material
- Comfortable light earloop



Item #	Description	Qty/Pack	Barcode
OPH-415	Kids 3D Disposable Face Mask (Medium)	20Pcs	 7 261894 304152

KIDS DISPOSABLE FACE MASK

3 protective Layers

- Sized for kids, making it comfy and fun!
- Light breathable material
- Comfortable light earloop

- Anti dust & fluid
- Comfortable & breathable
- Adjustable nose bridge
- Without glass fibers
- Hypoallergenic
- Very low resistance to breathing
- Nose bar adaptable
- High filtration capacity
- Perfect fitting



Item #	Description	Qty/Pack	Barcode
OPH-417	Kids Disposable Face Mask (Blue)	50Pcs	 7 261894 304176

KIDS DISPOSABLE FACE MASK

3 protective Layers

- Sized for kids, making it comfy and fun!
- Light breathable material
- Comfortable light earloop

- Anti dust & fluid
- Comfortable & breathable
- Adjustable nose bridge
- Without glass fibers
- Hypoallergenic
- Very low resistance to breathing
- Nose bar adaptable
- High filtration capacity
- Perfect fitting



Item #	Description	Qty/Pack	Barcode
OPH-418	Kids Disposable Face Mask (Pink)	50Pcs	 7 261894 304183

SPRAY HAND SANITIZER

Kills germs without water or towels



4 Different Scents:
Poppy, Hibiscus, Watermelon, & Cucumber Mint scented

70% Alcohol



Item #	Description	Qty/ml	Barcode
OPH-212	Spray Hand sanitizer	34ml	7 261894 302127



SPRAY HAND SANITIZER

Kills germs without water or towels

4 Different Scents:

Poppy, Hibiscus, Watermelon, & Cucumber Mint scented

70% Alcohol



Item #	Description	Qty/ml	Barcode
OPH-213	Spray Hand sanitizer	60ml	7 261894 302134

ANTIBACTERIAL HAND SANITIZER

Kills germs without water or towels

62% Alcohol



OPH - 203



OPH - 204



OPH - 209



OPH - 210



Item #	Description	Qty/ml	Barcode
OPH-203	Hand Sanitizer	53 ml	7 261894 302035
OPH-204	Hand Sanitizer	90 ml	7 261894 302042

Item #	Description	Qty/ml	Barcode
OPH-209	Hand sanitizer holder	1000 ml	7 261894 302097
OPH-210	Hand Sanitizer	1000 ml	7 261894 302103

ANTIBACTERIAL HAND SANITIZER

Kills germs without water or towels

75% Alcohol



OPH - 201




OPH - 202



OPH - 211



Item #	Description	Qty/ml	Barcode
OPH-201	Hand Sanitizer (75% Alcohol)	500 ml	7 261894 302011
OPH-202	Hand Sanitizer (75% Alcohol)	236 ml	7 261894 302028
OPH-211	Hand Sanitizer (75% Alcohol)	60 ml	7 261894 302110



62% Alcohol



OPH - 205 Cherry Blossom | Jasmine | Tropical Forest | Wood Out



OPH - 206 Cherry Blossom | Jasmine | Tropical Forest



- Cherry Blossom
- Jasmine
- Tropical Forest
- Strawberry
- Lemon
- Ocean

OPH - 207

Item #	Description	Qty/ml	Barcode
OPH-205	Hand sanitizer	29 ml	7 ¹¹ 261894 ¹¹ 302059
OPH-206	Hand sanitizer with silicone holder	29 ml	7 ¹¹ 261894 ¹¹ 302066
OPH-207	Hand sanitizer with cartoon silicone holder	35 ml	7 ¹¹ 261894 ¹¹ 302073

RUBBER HOT WATER BAGS

Made of 100 % rubber



OPH - 302



OPH - 299



OPH - 301



OPH - 303

Item #	Description	Qty	Barcode
OPH-299	Rubber Hot Water Bag	1 Pc	7 261894 302998
OPH-301	Rubber Hot Water Bag	1 Pc	7 261895 103013
OPH-302	Rubber Hot Water Bag	1 Pc	7 261895 103020
OPH-303	Rubber Hot Water Bag	1 Pc	7 261895 103037

2 liters
of water

PREMIUM QUALITY GAUZE BANDAGE

Non-Sterile Roll - 1 Dozen

- 40's/40's, 30 x 24
- Do not use if package is opened or damaged
- For single use only
- Keep in a cool and dry place



PREMIUM QUALITY GAUZE SWAB

Non-Sterile - 12 Ply

- 40's/40's, 26 x 18
- Do not use if package is opened or damaged
- For single use only
- Keep in a cool and dry place



Item #	Description	Size	Barcode
OPH-121	Non-Sterile Gauze Bandage	5cmx5yards	7 261894 301212
OPH-122	Non-Sterile Gauze Bandage	7.5cmx5yards	7 261894 301229
OPH-123	Non-Sterile Gauze Bandage	10cmx5yards	7 261894 301236
OPH-124	Non-Sterile Gauze Bandage	12cmx5yards	7 261894 301243

Item #	Description	Size	Barcode
OPH-155	Non-Sterile Gauze Swab	5 x 5	7 261894 301557
OPH-157	Non-Sterile Gauze Swab	7.5 x 7.5	7 261894 301571
OPH-160	Non-Sterile Gauze Swab	10 x 10	7 261894 301601

PREMIUM QUALITY STERILE GAUZE SWAB

Wound Dressing

- 40's/40's, 26 x 18
- Do not use if package is opened or damaged
- For single use only
- Keep in a cool and dry place

Sterile Paraffin Gauze BP



Paraffin gauze can be used for the treatment of burns, ulcers, skin grafts, and a variety of traumatic injuries.



Item #	Description	Size	Barcode
OPH-005	Sterile Gauze Swab	20 x 20 (10/box)	7 261894 300109
OPH-007	Sterile Gauze Swab	30 x 30 (10/box)	7 261894 300123
OPH-010	Sterile Gauze Swab	40 x 40 (10/box)	7 261894 300147
OPH-011	Sterile Paraffin Gauze BP	10 x 10 (10/box)	7 261894 200119

Item #	Description	Size	Barcode
OPH-105	Sterile Gauze Swab	20 x 20 (100/box)	7 261894 300116
OPH-107	Sterile Gauze Swab	30 x 30 (100/box)	7 261894 300079
OPH-110	Sterile Gauze Swab	40 x 40 (100/box)	7 261894 301106

SYRINGES FOR SINGLE USE

With Needles | Luer Lock

- STERILE
- Single use only
- Non-toxic
- Non-pyrogenic
- Gauge: 23G x 1 1/4" (0.6X32mm)
- Do not use if package is damaged

3CC/ML



Item #	Description	Qty/Pack	Barcode
OPH-125	Syringes for single use with needle (3cc)	100 Pcs	

SYRINGES FOR SINGLE USE

With Needles | Luer Lock

- STERILE
- Single use only
- Non-toxic
- Non-pyrogenic
- Gauge: 23G x 1 1/4" (0.6X32mm)
- Do not use if package is damaged

5CC/ML



Item #	Description	Qty/Pack	Barcode
OPH-126	Syringes for single use with needle (5cc)	100 Pcs	

INSULIN PEN NEEDLES

0.25 x 4mm | 31G



0.25 x 6mm | 31G



0.25 x 8mm | 31G



- Disposable, single use only
- Non-toxic, non-pyrogenic
- Fits most insulin delivery systems
- Secure attachment
- Smooth injection

Item #	Description	Qty/Pack	Barcode
OPH-130	Insulin Pen Needles 0.25 X 4mm 31g	100 Pcs	7 261894 301304
OPH-131	Insulin Pen Needles 0.25 X 6mm 31g	100 Pcs	7 261894 301311
OPH-132	Insulin Pen Needles 0.25 X 8mm 31g	100 Pcs	7 261894 301328

STERILE INSULIN SYRINGE FOR SINGLE USE

0.5CC/ML | 30 G



- STERILE
- Single use only
- Non-toxic
- Pyrogen free
- Gauge: 30g x 5/16" (8mm)
- Do not use if package is damaged



Item #	Description	Qty/Pack	Barcode
OPH-127	Sterile Insulin Syringe (0.5CC/ML)	100 Pcs	 7 261894 301274

1CC/ML | 30 G



- STERILE
- Single use only
- Non-toxic
- Pyrogen free
- Gauge: 30g x 5/16" (8mm)
- Do not use if package is damaged



Item #	Description	Qty/Pack	Barcode
OPH-128	Sterile Insulin Syringe (1CC/ML)	100 Pcs	 7 261894 301281



ALOE VERA PLASTERS WITH A HEALING ACTION

Antiseptic, elastic gauze with ALOE VERA extract; a well-known natural emollient, soothing and moisturizing agent. It protects the wound from external factors and enables natural healing.



Soothing, Moisturizing, Water-resistant, Breathable, Protective, Hypoallergenic, PVC and Latex free, Dermatologically tested

OPH - 425

30 Assorted Plasters



ECO-FRIENDLY PLASTERS WITH HIGH-THICKNESS GAUZE

A new line made with compostable support derived from sugar cane. Water resistant and breathable, with low environmental impact. Certified FSC. Individually wrapped and dermatologically tested.



100% pure cellulose gauze, Microperforated, Water-resistant, Breathable, Protective, Hypoallergenic, PVC and Latex free, Dermatologically tested

OPH - 436

30 Assorted Plasters



OPH - 426



OPH - 427



OPH - 428



OPH - 429

Item #	Description	Qty	Barcode
OPH-425	Aloe Vera plasters (assorted)	30	7 261894 304251
OPH-436	ECO FRIENDLY plasters (assorted)	30	7 261894 304367
OPH-426	PE Classic plasters medium size (7x2cm)	10	7 261894 304268
OPH-427	PE Classic plasters medium size (7x2cm)	20	7 261894 304275
OPH-428	PE Classic plasters (assorted)	30	7 261894 304282
OPH-429	PE Classic plasters medium size (7x2cm)	100	7 261894 304299

BANDAGES



- Corn cure plasters; with salicylic acid
- Relieves and removes corns
- Immediate pressure relief
- Quick and easy application
- Adhesive plaster holds pad in place
- Soft and protective
- Dermatologically tested

OPH - 422



OPH - 320



OPH - 324



OPH - 325



OPH - 326



OPH - 327



OPH - 328



OPH - 329

Item #	Description	Qty	Barcode
OPH-422	Corn cure plasters (with salicylic acid)	6	7 261 894 304220
OPH-320	Embossed clear-lite (Regular)	30	7 261 894 303209
OPH-324	Grid (Regular)	30	7 261 894 303247
OPH-325	Grid (Large +)	30	7 261 894 303254
OPH-326	Grid (Assorted)	36	7 261 894 303261
OPH-327	Complexion (Regular)	30	7 261 894 303278
OPH-328	Complexion (Large +)	30	7 261 894 303285
OPH-329	Complexion (Assorted)	36	7 261 894 303292

BANDAGES



OPH - 339



OPH - 343



OPH - 435



- 1 - Stick adhesive:
Long-lasting, unique adhesive helps bandage stay on
- 2 - Non-adhesive pad
Highly absorbent cushion protect cuts and scrapes without sticking to them
- 3 - Wrapper
Glazed paper, easy to open

Item #	Description	Qty	Barcode
OPH-339	Sheer Round (Regular)	50	7 261894 303391
OPH-343	Fabric Classic (Regular)	100	7 261894 303438
OPH-435	Aquastop waterproof plasters - XL (7.5x5cm)	5	7 261894 304350

WATERPROOF DRESSING

10 Pieces

- Infection protection
- Flexible and conformable
- Keeps the wound clean
- Antiseptic
- Hypoallergenic
- Sterile



Item #	Description	Qty/Pack	Barcode	Item #	Description	Qty/Pack	Barcode
OPH-358	Waterproof Dressing (5cm X 7.2cm)	10 Pcs	7 261894 303582	OPH-362	Waterproof Dressing (9cm X 15cm)	10 Pcs	7 261894 303629
OPH-359	Waterproof Dressing (8cm X 10cm)	10 Pcs	7 261894 303599	OPH-363	Waterproof Dressing (10cm X 20cm)	10 Pcs	7 261894 303636
OPH-360	Waterproof Dressing (8cm X 15cm)	10 Pcs	7 261894 303605	OPH-364	Waterproof Dressing (10cm X 25cm)	10 Pcs	7 261894 303643
OPH-361	Waterproof Dressing (9cm X 10cm)	10 Pcs	7 261894 303612	OPH-365	Waterproof Dressing (10cm X 35cm)	10 Pcs	7 261894 303650

WATERPROOF DRESSING

50 Pieces

- Infection protection
- Flexible and conformable
- Keeps the wound clean
- Antiseptic
- Hypoallergenic
- Sterile



OPH - 464



OPH - 465



OPH - 466



OPH - 467



OPH - 468



OPH - 469

Item #	Description	Qty/Pack	Barcode	Item #	Description	Qty/Pack	Barcode
OPH-464	Waterproof Dressing (5cm X 7.2cm)	50 Pcs	7 261894 304640	OPH-467	Waterproof Dressing (10cm X 20cm)	50 Pcs	7 261894 304671
OPH-465	Waterproof Dressing (8cm X 10cm)	50 Pcs	7 261894 304657	OPH-468	Waterproof Dressing (10cm X 25cm)	50 Pcs	7 261894 304688
OPH-466	Waterproof Dressing (9cm X 15cm)	50 Pcs	7 261894 304664	OPH-469	Waterproof Dressing (10cm X 35cm)	50 Pcs	7 261894 304695

NON-WOVEN DRESSING



OPH - 430

- White color
- Non-Woven
- Breathable
- For sensitive skin
- High-protection
- Hypoallergenic
- Soft and protective
- PVC and Latex free
- Dermatologically tested

Extra Sticky

Long-lasting, unique adhesive helps bandage stay on.

Non-Adhesive Pad

Highly absorbent cushion, protects cuts and scrapes without leaving an adhesive on the skin when removed



OPH - 433



OPH - 434



OPH - 431



OPH - 432

Item #	Description	Qty/Pack	Barcode
OPH-430	Non-Woven plasters (assorted)	30 Pc	7 261894 304305
OPH-431	Non-Woven plasters (XL size) (7.5x5cm)	5 Pc	7 261894 304312
OPH-432	Non-Woven plasters (XXXL size) (10x15cm)	6 Pc	7 261894 304329
OPH-433	STERILE Non-Woven plasters (XXXL size) (10x20cm)	6 Pc	7 261894 304336
OPH-434	STERILE Non-Woven plasters (XXXL size) (10x25cm)	6 Pc	7 261894 304343

NON-WOVEN DRESSING

10 Pieces

- Infection protection
- Flexible and conformable
- Keeps the wound clean
- Antiseptic
- Hypoallergenic
- Sterile



Non-Woven Dressing With Plastic Ring



Item #	Description	Qty/Pack	Barcode
OPH-373	Non Woven Dressing (10 X 20 cm)	10 Pcs	7 261894 303735
OPH-374	Non Woven Dressing (10 X 25 cm)	10 Pcs	7 261894 303742
OPH-375	Non Woven Dressing (10 X 35 cm)	10 Pcs	7 261894 303759
OPH-376	Non Woven Dressing (8 X 10 cm)	10 Pcs	7 261894 303766

Item #	Description	Qty/Pack	Barcode
OPH-377	Non Woven Dressing (9 X 15 cm)	10 Pcs	7 261894 303773
OPH-412	Non Woven Dressing (5 X 7.2 cm)	10 Pcs	7 261894 304121
OPH-421	Non Woven Dressing with plastic ring (10 X 10 cm)	10 Pcs	7 261894 304213

NON-WOVEN DRESSING

50 Pieces

- Infection protection
- Flexible and conformable
- Keeps the wound clean
- Antiseptic
- Hypoallergenic
- Sterile



OPH - 470



OPH - 471



OPH - 472



OPH - 473



OPH - 474

Item #	Description	Qty/Pack	Barcode
OPH-470	Non-Woven Dressing (10cm X 20cm)	50 Pcs	7 261894 304701
OPH-471	Non-Woven Dressing (10cm X 25cm)	50 Pcs	7 261894 304718
OPH-472	Non-Woven Dressing (10cm X 35cm)	50 Pcs	7 261894 304725

Item #	Description	Qty/Pack	Barcode
OPH-473	Non-Woven Dressing (8cm X 10cm)	50 Pcs	7 261894 304732
OPH-474	Non-Woven Dressing (9cm X 15cm)	50 Pcs	7 261894 304749

WOUND CLOSURE

Adhesive strips

- Comfortable adhesive strips used to secure, close and support small cuts & wounds
- They stay firmly in place, allowing the skin to breathe as the wound heals
- Use wound closure strips as directed by your physician to treat small cuts and wounds
- Wound closure strips may also be used to support skin sutures and staples or after their removal, for wound support



EYE PATCH

Sterile With Adhesive Edge

- Soft
- For use in the treatment of eye injuries
- Post to the application of ointments
- Post operative eye care
- Non-stick patch won't stick to wounds
- Unique, long-lasting adhesive edge ensures that the eye patch stays in place and seal out dirt and germs

For Kids



Item #	Description	Qty/Pack	Barcode
OPH-345	Wound Closure (6.35 x 101 mm)	30 Pcs	7 261894 303452
OPH-369	Wound Closure (12 x 101 mm)	15 Pcs	7 261894 303698
OPH-423	Wound Closure (12.7 x 108 mm)	18 Pcs	7 261894 304237

Item #	Description	Qty/Pack	Barcode
OPH-346	Eye patch for kids (63 x 57 mm)	10 Pcs	7 261894 303469
OPH-370	Eye patch (85 x 60 mm)	10 Pcs	7 261894 303704

FIRST AID KIT

Used by Healthcare professionals



- 6 BZK wipes (1.97"x1.97")
- 4 Alcohol cleansing pads (1.97"x1.97")
- 25 2.99" x 0.75" PE bandage - Skin color
- 1 2.99" x 0.98" PE bandage - Skin color
- 20 1.57" x 0.39" PE bandage - Skin color
- 3 1.50" x 1.50" PE bandage - Skin color
- 1 3.94" x 1.97" PE bandage - Skin color
- 2 1.97" x 1.50" Flexible fabric bandage
- 2 2.01" x 1.77" Flexible fabric bandage
- 1 2.76" x 3.94" Non-Woven Dressing
- 1 3.35" x 2.28" eye pad

- 1 6" x 3/4" Finger splints
- 1 2" x 2" Gauze pad
- 1 3" x 3" Gauze pad
- 1 0.49" x 180" Latex free tape
- 2 1.97" x 157.48" Sterile gauze roll
- 12 3" Cotton tipp - Single Stick
- 1 4.53" Scissors
- 1 4.92" Tweezers
- 6 Pins - Middle size
- 2 Large gloves (in a ziplock)
- 1 Plastic case

OPH - 378

95 pieces



- 4 BZK wipes (1.97"x1.97")
- 4 Alcohol cleansing pads (1.97"x1.97")
- 25 2.99" x 0.75" PE bandage - Regular
- 2 2.99" x 0.98" PE bandage - Large
- 30 1.57" x 0.39" PE bandage - Small
- 2 1.50" x 1.50" PE bandage - Patch
- 1 3.94" x 1.97" PE bandage - Extra Large
- 2 1.97" x 1.50" Knuckle fabric bandage
- 2 2.01" x 1.77" Fingerstrip fabric bandage
- 1 2.76" x 0.47" Butterfly Closure
- 1 6" x 3/4" Finger splints

- 1 3.35" x 2.28" Sterile eye pad
- 1 2" x 2" Gauze pads
- 1 3" x 3" Gauze pads
- 1 5" x 9" Sterile Trauma pad
- 1 0.49" x 180" Latex free tape
- 1 1.97" x 157.48" Sterile gauze roll
- 12 3" Cotton tipp - Single Stick
- 2 Large PVC gloves
- 1 4.53" Scissors
- 1 4.92" Tweezers
- 1 Carry anywhere plastic case
- 1 First aid guide

OPH - 379

98 pieces



Item #	Description	Qty/Pack	Barcode
OPH-378	First Aid Kit	95 Pcs	
OPH-379	First Aid Kit	98 Pcs	



- Treats cuts, scrapes, bruises, minor aches and pains
- Family first aid kit contains an ideal mix of first aid items to look after your family at home and on the go.

KIDS BANDAGES



OPH - 386

Luminous



OPH - 437



OPH - 438



OPH - 383



OPH - 384



OPH - 385



OPH - 387

Item #	Description	Qty/Pack	Barcode
OPH-386	Kids Bandages (Luminous)	20 Pcs	7 261894 303865
OPH-437	Kids Bandages (Ocean)	20 Pcs	7 261894 304374
OPH-438	Kids Bandages (Emoji)	20 Pcs	7 261894 304381
OPH-383	Kids Bandages	16 Pcs	7 261894 303834
OPH-384	Kids Bandages	16 Pcs	7 261894 303841
OPH-385	Kids Bandages	16 Pcs	7 261894 303858
OPH-387	Kids Bandages	12 Pcs	7 261894 303872

BANDAGE ROLL | LATEX FREE

The bandage roll secures the bandage to the wound by sticking to itself, not to skin or hair, and won't leave a sticky residue when changing the bandage. Breathable and fully flexible

USE SELF-ADHERING WRAP TO:

- PROTECT skin
- SECURE dressings, catheters
- SUPPORT sprains/strains, splints, compresses for edema



OPH - 461
(5 CM X 4.5 M)

Item #	Description	Qty/Pack	Barcode - Blue	Barcode - Green
OPH-461	Bandage Roll (5cm X 4.5m)	24 Pc	 7 261893 304610	 7 261890 304613
OPH-462	Bandage Roll (7.5cm X 4.5m)	16 Pc	 7 261893 304627	 7 261890 304620

Medical Grade Quality

Custom Weave
No Scissors Needed

Proprietary Adhesive
Does Not Pull Out Your Hair



Natural latex



OPH - 462
(7.5 CM X 4.5 M)

BANDAGE ROLL | LATEX FREE

The bandage roll secures the bandage to the wound by sticking to itself, not to skin or hair, and won't leave a sticky residue when changing the bandage. Breathable and fully flexible

USE SELF-ADHERING WRAP TO:

- PROTECT skin
- SECURE dressings, catheters
- SUPPORT sprains/strains, splints, compresses for edema



OPH - 459
(7.5 CM X 4.5 M)



OPH - 458
(5 CM X 4.5 M)



OPH - 460
(10 CM X 4.5 M)

Item #	Description	Qty/Pack	Barcode - Beige	Barcode - White
OPH-458	Bandage Roll (5cm X 4.5m)	24 Pc	 7 261899 304584	 7 261892 304581
OPH-459	Bandage Roll (7.5cm X 4.5m)	16 Pc	 7 261899 304591	 7 261892 304598
OPH-460	Bandage Roll (10cm X 4.5m)	12 Pc	 7 261899 304607	 7 261892 304604

TRANSPARENT PE TAPE



OPH - 451 | 1.25cm x 9.14m



OPH - 452 | 2.5cm x 9.14m



OPH - 453 | 5cm x 9.14m



- Ultra-thin, super soft and comfortable to wear
- Waterproof, strong glue & breathable



A transparent tape with a polyethylene backing and an acrylic adhesive that helps make repairs less visible. The polyethylene backing is conformable and provides good chemical and solvent resistance

Item #	Description	Qty/Pack	Barcode
OPH-451	Transparent PE tape (1.25cm x 9.14m)	24 Rolls	7 261894 304510
OPH-452	Transparent PE tape (2.5cm x 9.14m)	12 Rolls	7 261894 304527
OPH-453	Transparent PE tape (5cm x 9.14m)	6 Rolls	7 261894 304534

NON STERILE ELASTIC BANDAGES



OPH - 391
5cm x 4m



OPH - 392
7.5cm x 4m



OPH - 393
10cm x 4m



OPH - 403
36inch x 50yds



OPH - 394
5cm x 4m



OPH - 395
7.5cm x 4m



OPH - 396
10cm x 4m



OPH - 455
5cm x 4.5m



OPH - 456
7.5cm x 4.5m



OPH - 457
10cm x 4.5m

Item #	Description	Qty/Pack	Barcode
OPH-391	Non Sterile Elastic Bandage (5cm x 4m)	1 Pc	7 261894 303919
OPH-392	Non Sterile Elastic Bandage (7.5cm x 4m)	1 Pc	7 261894 303926
OPH-393	Non Sterile Elastic Bandage (10cm x 4m)	1 Pc	7 261894 303933
OPH-403	Absorbent Gauze Roll (36inch x 50yds)	1 Pc	7 261894 304039
OPH-394	Non Sterile Elastic Bandage (5cm x 4m)	1 Pc	7 261894 303940

Item #	Description	Qty/Pack	Barcode
OPH-395	Non Sterile Elastic Bandage (7.5cm x 4m)	1 Pc	7 261894 303957
OPH-396	Non Sterile Elastic Bandage (10cm x 4m)	1 Pc	7 261894 303964
OPH-455	Non Sterile Elastic Bandage (5cm x 4.5m)	1 Pc	7 261894 304558
OPH-456	Non Sterile Elastic Bandage (7.5cm x 4.5m)	1 Pc	7 261894 304565
OPH-457	Non Sterile Elastic Bandage (10cm x 4.5m)	1 Pc	7 261894 304572





OPH - 397
5cm x 4.5m



OPH - 398
7.5cm x 4.5m



OPH - 399
10cm x 4.5m



OPH - 400
5cm x 4.5m



OPH - 401
7.5cm x 4.5m



OPH - 402
10cm x 4.5m



OPH - 390
10cm x 4.5m



OPH - 389
7.5cm x 4.5m



OPH - 388
5cm x 4.5m



Item #	Description	Qty/Pack	Barcode
OPH-397	High Elastic Bandage (5cm x 4.5m)	1 Pc	7 261894 303971
OPH-398	High Elastic Bandage (7.5cm x 4.5m)	1 Pc	7 261894 303988
OPH-399	High Elastic Bandage (10cm x 4.5m)	1 Pc	7 261894 303995
OPH-400	High Elastic Bandage (5cm x 4.5m)	1 Pc	7 261894 304008
OPH-401	High Elastic Bandage (7.5cm x 4.5m)	1 Pc	7 261894 304015

Item #	Description	Qty/Pack	Barcode
OPH-402	High Elastic Bandage (10cm x 4.5m)	1 Pc	7 261894 304022
OPH-388	Non Sterile Elastic Bandage (5cm x 4.5m)	1 Pc	7 261894 303889
OPH-389	Non Sterile Elastic Bandage (7.5cm x 4.5m)	1 Pc	7 261894 303896
OPH-390	Non Sterile Elastic Bandage (10cm x 4.5m)	1 Pc	7 261894 303902

PLASTER OF PARIS BANDAGES



OPH - 404
7.5cm x 2.7 m



OPH - 405
10cm x 2.7 m



OPH - 406
15cm x 2.7 m

- For orthopedic fixation, deformity correction, inflammation limb brakes, osteoarticular surgery, limb fixation and mold making
- Provides excellent stability and conformability

Item #	Description	Qty/Pack	Barcode
OPH-404	Plaster of Paris Bandage (7.5cm X 2.7m)	1 Pc	 7 ²⁶¹⁸⁹⁴ 304046
OPH-405	Plaster of Paris Bandage (10cm X 2.7m)	1 Pc	 7 ²⁶¹⁸⁹⁴ 304053
OPH-406	Plaster of Paris Bandage (15cm X 2.7m)	1 Pc	 7 ²⁶¹⁸⁹⁴ 304060



NONWOVEN MEDICAL TAPE

- Ultra-thin, super soft and comfortable to wear | Waterproof & strong glue & breathable
- Use for minor wound such as cuts, scrapes, abrasions & blisters



OPH - 351 | 1.25cm x 9.14m



OPH - 352 | 2.5cm x 9.14m



OPH - 353 | 5cm x 9.14m



OPH - 354 | 7.5cm x 9.14m

Item #	Description	Qty/Pack	Barcode
OPH-351	Nonwoven Medical Tape (1.25cm X 9.14m)	24 Rolls	7 261894 303513
OPH-352	Nonwoven Medical Tape (2.5cm X 9.14m)	12 Rolls	7 261894 303520

Item #	Description	Qty/Pack	Barcode
OPH-353	Nonwoven Medical Tape (5cm X 9.14m)	06 Rolls	7 261894 303537
OPH-354	Nonwoven Medical Tape (7.5cm X 9.14m)	06 Rolls	7 261894 303544

SURGICAL COTTON

Ideal for all your baby care, skin care and family needs

- For external use only
- 100% Pure cotton
- Recommended for all medical and dermatological uses.



MAKE-UP REMOVER PADS

Ideal for removing make-up, nail polish, applying creams and lotions, cleaning eyes, ears, skin of babies and children.

Enjoy cleaning your face and eyes with OPTIMAL pads and feel the real softness of cotton.



20 FREE PADS

Item #	Description	Qty	Barcode
OPH-074	Surgical Cotton	50 g.	7 261 895 100074
OPH-081	Surgical Cotton	100 g.	7 261 895 100081
OPH-098	Make-up Remover pads (round)	80 Pcs	7 261 895 100098
OPH-096	Make-up Remover pads (round)	120 Pcs	7 261 895 100968

WOOL BALLS

Soft, Gentle, Hygienic



- Ideal for all your baby care, skin care and family needs
- For external use only
- 100% Pure cotton
- Recommended for all medical and dermatological uses.

Item #	Description	Qty	Barcode
OPH-111	White Cotton Wool balls	100	

COTTON BUDS

Recommended for new borns



OPH - 055
55 pcs

Recommended for new borns



OPH - 056
80 pcs

Recommended for new borns



OPH - 057
108 pcs

Item #	Description	Qty	Barcode
OPH-055	Cotton Buds	55pcs	
OPH-056	Cotton Buds	80pcs	
OPH-057	Cotton Buds	108pcs	

COTTON BUDS

- For external use only
- 100% Pure cotton
- Recommended for all medical and dermatological uses.



OPH - 100
100 pcs



OPH - 200
200 pcs



OPH - 300
300 pcs

Item #	Description	Qty	Barcode
OPH-100	Cotton Buds	100pcs	7 261895 101002
OPH-200	Cotton Buds	200pcs	7 261895 102009
OPH-300	Cotton Buds	300pcs	7 261895 103006
OPH-290	Wood Cotton Buds	300pcs	7 261894 302905

Wood Cotton Buds


OPH - 290
300 pcs



OVULATION TEST

Over 99.9% Accurate | predicts ovulation 24 - 48 hours in advance results in 5 minutes


- The result window will show a clear result. Pregnatest is over 99% accurate.
- Hold the color change tip in the urine stream for 10 secs.
- In 5 minutes the line in the end window will show that the test has been completed and your result is ready. The result window will show a clear result.





PREGNANCY TESTS

Fast, Clear and Accurate

One Step HCG Urine Test



- Use the pipette to withdraw urine specimen, dispense 4 drops into the sample well.
- Read result within 5 minutes.



- This device is for self-test.
- Store at 4~30°C.
- Keep out of reach of children.
- Pregnatest is a trademark.

- The result window will show a clear result. Pregnatest is over 99% accurate.
- Hold the color change tip in the urine stream for 10 secs
- At 3 minutes the line in the end window will show that the test has been completed and your result is ready. The result window will show a clear result

Item #	Description	Qty	Barcode
OPM-144	Ovulation test	5	7 261895 1101446

Item #	Description	Qty	Barcode
OPM-142	One Step HCG Urine Test	1	7 261895 1100142
OPM-143	Pregnancy test	1	7 261895 1101439

DISPOSABLE UNDERPAD

Waterproof backing Hygienic

- Soft, Breathable, and Absorbent
- Highly absorbent
- Non-woven very soft top sheet
- Waterproof back sheet 100% Polyethylene
- Prevention of fluid leakage to help keep surfaces dry
- Durable



LARGE (60cm x 90cm)



For Men and Women



✓ **Normal**



Item #	Description	Qty/Pack	Barcode
OPH-501	Disposable Underpad	15Pcs	7 261894 305012

FEMININE INCONTINENCE PAD

Suprema-Gyno

For Heavy Flow: Breathable | Discreet | Absorbent

- Latex Free
- Reduction of unpleasant smell
- Super Absorbent Core
- Skin breathability
- Strong adhesive tape
- Easy and reliable fixing to underwear



For Women



Super absorbent Polymer



Maximum Absorbency

Item #	Description	Qty/Pack	Barcode
OPH-500	Feminine Incontinence Pad	14Pcs	7 261894 305005

REFRESHING WIPES

With The Soothing Touch Of Aloe Vera



0%
Alcohol
Paraben
Sulphate
Chlorine
Phenoxyethanol

- Soft and Gentle
- Hypoallergenic
- Enriched with Vitamin E and Aloe Vera



72
Scented
Wipes

Item #	Description	Qty/Pack	Barcode
OPH-060	Refreshing wipes	72 wipes	 7 261895 100609

REFRESHING WIPES

With The Soothing Touch Of Aloe Vera



0%
Alcohol
Paraben
Sulphate
Chlorine
Phenoxyethanol

- Soft and Gentle
- Hypoallergenic
- Enriched with Vitamin E and Aloe Vera



52
Scented
Wipes

Item #	Description	Qty/Pack	Barcode
OPH-061	Refreshing Back To School wipes	52 wipes	 7 261895 100616

ANTIBACTERIAL WIPES

Soft and Fresh



15 wipes



72 wipes

- Alcohol free
- With Aloe Vera/ Chamomile extract.
- Microbiologically tested.
- Keeps natural Ph-balance of your skin

Item #	Description	Qty/Pack	Barcode
OPH-015	Antibacterial wipes	15 Pcs	7 261894 300154
OPH-065	Antibacterial wipes	72 Pcs	7 261895 100654

HAND AND BODY WIPES

Soft and Fresh



- Cleanses & refreshes in just one step: Wipes for hands and body.
- Perfect to go! Dermatologically tested.
- ADVICE: do not apply on irritated or wounded skin. Keep away from children. Reseal sachet after use.

Item #	Description	Qty/Pack	Barcode
OPH-062	Hand and body wipes	15 wipes	7 261895 1100623

HYDRATING HAND CREAM

Moisturizing hand cream with Aloe Vera, extract of Asiatic Centella and Hamamelis. Instantly absorbs into skin, moisturizes and protects it from the harmful irritants.



SOFTENING FOOT CREAM

Softening foot cream with Shea Butter and Vitamin E, softens and moisturizes skin.



Item #	Description	Qty/ml	Barcode
OPH-1600	Hydrating Hand Cream	125ml	7 261894 316001
OPH-1601	Softening Foot Cream	125ml	7 261894 316018



- Paraben Free
- Dermatologically tested



SUNSCREEN PRODUCTS

OPTIMAL® Sunscreen locks active ingredients in place. It decreases the risk of skin cancer and early skin aging caused by the sun, and combines a highly-effective UVA/UVB system which is absorbed very quickly, without leaving any perceivable residue on the skin.

Water resistant | Paraben free | UVA UVB protection



SPRAY SUNSCREEN
(SPF 50+) - 250 ml



SUNSCREEN
(SPF 50+) - 250 ml



SUNSCREEN
(SPF 50+) - 125 ml

Item #	Description	Qty/ml	Barcode
OPH-1606	Spray Sunscreen SPF 50+	250ml	7 261894 316063
OPH-1603	Sunscreen SPF 50+	250ml	7 261894 316032
OPH-1605	Sunscreen SPF 50+	125ml	7 261894 316056

COCONUT SUN OIL WITH VITAMIN E (SPF 6)

Optimal tanning oil SPF 6 is a transparent sun oil that protects your skin against UVA and UVB rays. Water resistant. The formula with Vitamin E leaves a smooth and silky skin feeling.



Item #	Description	Qty/ml	Barcode
OPH-1604	Coconut Sun oil SPF 6	200ml	7 261894 316049



OPTIMAL Medics s.a.r.l.

Health Care Products

Distributed by:

UNIMARC S.A.R.L.

Head Office:

Bd des Eplatures 64, 2300 La Chaux-de-Fonds, Switzerland

Tel: +41 762 906 320, E-mail: info@unimarc.ch

Website: www.unimarc.ch

OPTIMAL MEDICS S.A.R.L.

Head Office:

Horizon Building, Doha Entrance, Choueifat, Lebanon

Mob: +961 70 555 445, Tel: +961 5 812 850

E-mail: info@optimalmedics.com, Website: www.optimalmedics.com



Visit us on:



Optimal



Optimal_medics